

anti- AFAP1L2 antibody

Product Information

Catalog No.: FNab00194

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

May play a role in a signaling cascade by enhancing the kinase activity of SRC.Contributes to SRC-regulated transcription activation.

Immunogen information

Immunogen: actin filament associated protein 1-like 2

Synonyms: AFAP1 like protein 2, AFAP1L2, FLJ14564, KIAA1914, XB130

Observed MW: 110 kDa UniprotID: Q8N4X5

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB

Recommended dilution: WB: 1:200-1:1000; IP: 1:200-1:2000

Image:


1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





mouse thymus tissue were subjected to SDS PAGE followed by western blot with FNab00194(AFAP1L2 Antibody) at dilution of 1:300