

anti- AGBL5/CCP5 antibody

Product Information

Catalog No.:	FNab00208
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	ATP/GTP binding protein-like 5
Synonyms:	AGBL5, ATP/GTP binding protein like 5, CCP5, FLJ21839
Observed MW:	75 kDa
Uniprot ID :	Q8NDL9

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Fusion protein expressed in Baculovirus were subjected to SDS PAGE followed by western blot with FNab00208(AGBL5 Antibody) at dilution of 1:500