

anti- AGPAT3 antibody

Product Information

Catalog No.: FNab00216

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Converts lysophosphatidic acid(LPA) into phosphatidic acid by incorporating an acyl moiety at the sn-2 position of the glycerol backbone. Acts on LPA containing saturated or unsaturated fatty acids C16:0-C20:4 at the sn-1 position using C18:1, C20:4 or C18:2-CoA as the acyl donor. Also acts on lysophosphatidylcholine, lysophosphatidylinositol and lysophosphatidylserine using C18:1 or C20:4-CoA(PubMed:21173190). Has a preference for arachidonoyl-CoA as a donor. Has also a modest lysophosphatidylinositolacyltransferase(LPIAT) activity, converts lysophosphatidylinositol(LPI) into phosphatidylinositol(By similarity).

Immunogen information

Immunogen: 1-acylglycerol-3-phosphate O-acyltransferase 3

Synonyms: 1 AGP acyltransferase 3, 1 AGPAT 3, AGPAT3, LPAAT gamma,

LPAAT GAMMA1

Observed MW: 36-40 kDa UniprotID: Q9NRZ7

Application

1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com




Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:5000

Image:



mouse testis tissue were subjected to SDS PAGE followed by western blot with

FNab00216(AGPAT3 Antibody) at dilution of

1:1000