

anti- AGPAT6 antibody

Product Information

Catalog No.:	FNab00217
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Background	

Duckgrou

None

Immunogen information

Immunogen:	1-acylglycerol-3-phosphate O-acyltransferase 6(lysophosphatidic acid acyltransferase, zeta)
Synonyms:	1 AGP acyltransferase 6, 1 AGPAT 6, AGPAT6, DKFZp586M1819, GPAT4, LPAAT zeta, LPAATZ, TSARG7
Observed MW:	48-50 kDa
UniprotID :	Q86UL3


Application

Reactivity:	Human, Mouse, Rat	
Tested Application:	ELISA, WB, IHC, IP	
Recommended dilution: WB: 1:500-1:5000; IP: 1:500-1:5000; IHC: 1:20-1:200		
Image:		

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan,Hubei, China(430206)Tel :(0086)027-87384275Fax: (0086)027-87800889 www.fn-test.com





Immunohistochemistry of paraffin-embedded human testis using FNab00217(AGPAT6 antibody) at dilution of 1:100


IP Result of anti-AGPAT6 (IP:FNab00217, 4ug; Detection:FNab00217 1:1000) with mouse brain tissue lysate 2640ug.

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com