

anti- AGTRAP antibody

Product Information

Catalog No.: FNab00223

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Appears to be a negative regulator of type-1 angiotensin II receptor-mediated signaling by regulating receptor internalisation as well as mechanism of receptor desensitization such as phosphorylation. Induces also a decrease in cell proliferation and angiotensin II-stimulated transcriptional activity.

Immunogen information

Immunogen: angiotensin II receptor-associated protein

Synonyms: AGTRAP, ATRAP

Observed MW: 18 kDa Uniprot ID: Q6RW13

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, IHC

Recommended dilution: IF: 1:20-1:200; IHC: 1:10-1:100

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>

1



Image:



Immunohistochemical of paraffin-embedded human pancreas cancer using FNab00223(AGTRAP antibody) at dilution of 1:10