

# anti- AGXT antibody

#### **Product Information**

Catalog No.: FNab00224

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

None

## Immunogen information

Immunogen: alanine-glyoxylate aminotransferase

Synonyms: AGT, AGT1, AGXT, AGXT1, PH1, SPAT, SPT, TLH6

Observed MW: 43 kDa UniprotID: P21549

## **Application**

Reactivity: Human, Mouse

Tested Application: ELISA, IHC, WB

Recommended dilution: WB: 1:500-1:5000; IHC: 1:20-1:200

Image:


#### 1

#### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Immunohistochemistry of paraffin-embedded human liver slide using FNab00224(AGXT Antibody) at dilution of 1:50



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab00224(AGXT antibody) at dilution of 1:1000