

anti- AHI1 antibody

Product Information

Catalog No.: FNab00229

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24

months (Avoid repeated freeze / thaw cycles.)

Background

Component of the tectonic-like complex, a complex localized at the transition zone of primary cilia and acting as a barrier that prevents diffusion of transmembrane proteins between the cilia and plasma membranes.

Immunogen information

Immunogen: Abelson helper integration site 1

Synonyms: AHI 1, AHI1, dJ71N10.1, DKFZp686J1653, FLJ14023, FLJ20069,

JBTS3, Jouberin, ORF1

Observed MW: 137 kDa Uniprot ID: Q8N157

Application

Reactivity: Human, Mouse Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:5000

Image:


1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





HeLa cells were subjected to SDS PAGE followed by western blot with FNab00229(AHI1 antibody) at dilution of 1:1000