

anti- AHSA1 antibody

Product Information

Catalog No.: FNab00232

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Cochaperone that stimulates HSP90 ATPase activity(By similarity). May affect a step in the endoplasmic reticulum to Golgi trafficking.

Immunogen information

Immunogen: AHA1, activator of heat shock 90kDa protein ATPase homolog 1(yeast)

Synonyms: AHA1, AHSA1, C14orf3

Observed MW: 38 kDa
UniprotID: 095433

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, IHC, IF, WB, IP

Recommended dilution: WB: 1:200-1:2000; IF: 1:20-1:200; IHC: 1:20-1:200; IP: 1:200-1:2000

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>






Immunohistochemistry of paraffin-embedded human gliomas using FNab00232(AHSA1 antibody) at dilution of 1:50



IP Result of anti-AHSA1 (IP:FNab00232, 4ug; Detection:FNab00232 1:500) with mouse brain tissue lysate 2640ug.



A549 cells were subjected to SDS PAGE followed by western blot with FNab00232(AHSA1 Antibody) at dilution of 1:600

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel: (0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>



Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel: (0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>