

## anti- AKT1 antibody

### Product Information

Catalog No.:	FNab00269
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	v-akt murine thymoma viral oncogene homolog 1
Synonyms:	AKT, AKT1, PKB, PKB ALPHA, PRKBA, Protein kinase B, Proto oncogene c Akt, RAC, RAC ALPHA, RAC PK alpha
Observed MW:	62 kDa
Uniprot ID :	P31749

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB, IF
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human breast cancer using FNab00269(AKT1 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00269(AKT1 antibody) at dilution of 1:500