

# anti- ALDH1A1 antibody

#### **Product Information**

Catalog No.: FNab00283

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

Binds free retinal and cellular retinol-binding protein-bound retinal. Can convert/oxidize retinaldehyde to retinoic acid(By similarity).

### **Immunogen information**

Immunogen: aldehyde dehydrogenase 1 family, member A1

Synonyms: ALDC, ALDH E1, ALDH1, ALDH11, ALDH1A1, ALHDII, PUMB1,

RALDH 1, RALDH1, Retinal dehydrogenase 1

Observed MW: 55 kDa UniprotID: P00352

# **Application**

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500-1:5000; IHC: 1:20-1:200

Image:


#### 1

#### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Immunohistochemistry of paraffin-embedded human liver using FNab00283(ALDH1A1 antibody) at dilution of 1:50



human liver tissue were subjected to SDS PAGE followed by western blot with FNab00283(ALDH1A1 antibody) at dilution of 1:1000