

anti- ALDH1B1 antibody

Product Information

Catalog No.: FNab00287

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

ALDHs play a major role in the detoxification of alcohol-derived acetaldehyde. They are involved in the metabolism of corticosteroids, biogenic amines, neurotransmitters, and lipid peroxidation.

Immunogen information

Immunogen: aldehyde dehydrogenase 1 family, member B1

Synonyms: Aldehyde dehydrogenase 5, ALDH1B1, ALDH5, ALDHX

Observed MW: 56 kDa Uniprot ID: P30837

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB, IHC, IP,IF

Recommended dilution: WB: 1:500-1:3000; IHC: 1:20-1:200; IF: 1:20-1:200; IP: 1:500-1:1000

Image:


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)


Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>

1





Immunohistochemistry of paraffin-embedded human testis using FNab00287(ALDH1B1 antibody) at dilution of 1:100




IP Result of anti-ALDH1B1 (IP: FNab00287, 4ug; Detection: FNab00287 1:500) with mouse brain tissue lysate 4000ug.

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





HepG2 cells were subjected to SDS PAGE followed by western blot with FNab00287(ALDH1B1 antibody) at dilution of 1:1000

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>