

# anti- ALDH1L2 antibody

#### **Product Information**

Catalog No.:	FNab00290
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

## Background

ALDH1L2(Aldehyde dehydrogenase family 1 member L2) is also named as mitochondrial 10-FTHFDH, mtFDH and belongs to the aldehyde dehydrogenase family and ALDH1L subfamily. The ALDH1L2 gene has three transcriptional variants with the molecular weight of 102 kDa, 42 kDa and 89 kDa(PMID:19823103).

#### Immunogen information

Immunogen:	aldehyde dehydrogenase 1 family, member L2
Synonyms:	ALDH1L2, DKFZp686A16126, DKFZp686M064, FLJ38508
Observed MW:	102 kDa and 89 kDa
UniprotID :	Q3SY69

## Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution	:WB: 1:200-1:2000; IHC: 1:20-1:200

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206) Tel :( 0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com



## Image:



-100kd

-70kd

Immunohistochemistry of paraffin-embedded human pancreas slide using FNab00290(ALDH1L2 Antibody) at dilution of 1:50

mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab00290(ALDH1L2 antibody) at dilution of 1:1000

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>