

anti- ALDH3B2 antibody

Product Information

Catalog No.: FNab00295

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Oxidizes medium and long chain aldehydes into non-toxic fatty acids.

Immunogen information

Immunogen: aldehyde dehydrogenase 3 family, member B2

Synonyms: ALDH3B2, Aldehyde dehydrogenase 8, ALDH8

Observed MW: 43 kDa UniprotID: P48448

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB, IHC, IF

Recommended dilution: WB: 1:500-1:5000; IHC: 1:20-1:200; IF: 1:10-1:100

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human adrenal gland using FNab00295(ALDH3B2 antibody) at dilution of 1:50



mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab00295(ALDH3B2 antibody) at dilution of 1:500