

anti- ALDH5A1 antibody

Product Information

Catalog No.: FNab00297

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Catalyzes one step in the degradation of the inhibitory neurotransmitter gamma-aminobutyric acid(GABA).

Immunogen information

Immunogen: aldehyde dehydrogenase 5 family, member A1

Synonyms: ALDH5A1, SSADH, SSDH
Observed MW: 57 kDa, observed 51 kDa

UniprotID: P51649

Application

Reactivity: Human, Mouse, Rat Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500-1:5000; IHC: 1:20-1:200

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Immunohistochemistry of paraffin-embedded human liver tissue slide using FNab00297(ALDH5A1 Antibody) at dilution of 1:50



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with FNab00297(ALDH5A1 antibody) at dilution of 1:500