

anti- ALDH6A1 antibody

Product Information

Catalog No.:	FNab00298
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Plays a role in valine and pyrimidine metabolism. Binds fatty acyl-CoA.

Immunogen information

Immunogen:	aldehyde dehydrogenase 6 family, member A1
Synonyms:	ALDH6A1, MMSADHA, MMSDH
Observed MW:	58 kDa
Uniprot ID :	Q02252

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:2000-1:10000; IHC: 1:50-1:500
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel :(0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com





Immunohistochemistry of paraffin-embedded human kidney using FNab00298(ALDH6A1 antibody) at dilution of 1:50



HEK293 cells were subjected to SDS PAGE followed by western blot with FNab00298(ALDH6A1 antibody) at dilution of 1:5000

2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)