

anti- ALDH7A1 antibody

Product Information

Catalog No.: FNab00299

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Multifunctional enzyme mediating important protective effects. Metabolizes betaine aldehyde to betaine, an important cellular osmolyte and methyl donor. Protects cells from oxidative stress by metabolizing a number of lipid peroxidation-derived aldehydes. Involved in lysine catabolism.

Immunogen information

Immunogen: aldehyde dehydrogenase 7 family, member A1

Synonyms: ALDH7A1, Alpha AASA dehydrogenase, Antiquitin 1, ATQ1, Betaine

aldehyde dehydrogenase, EPD, FLJ11738, FLJ92814, P6c

dehydrogenase, PDE

Observed MW: 55 kDa UniprotID: P49419

Application

Reactivity: Human, Mouse, Rat Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500-1:5000: IHC: 1:20-1:200

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com


1



Image:



Immunohistochemistry of paraffin-embedded human ovary tumor using FNab00299(ALDH7A1 antibody) at dilution of 1:50



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab00299(ALDH7A1 antibody) at dilution of 1:1000