

anti- Alpha galactosidase A antibody

Product Information

Catalog No.: FNab00330

Size: 100µ g Form: liquid

Purification: Protein A+G purification

Purity: ≥95% as determined by SDS-PAGE

Host: Mouse IsoType: IgG2a

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

GLA, also named as Melibiase, Agalsidase and Alpha-galactosidase A, belongs to the glycosyl hydrolase 27 family. It hydrolyzes terminal, non-reducing alpha-D-galactose residues in alpha-D-galactosides, including galactose oligosaccharides, galactomannans and galactolipids. Fabry disease is an X-linked lysosomal storage disorder resulting from the deficient activity of GLA. Enzyme replacement therapy (ERT) with GLA is currently the most effective therapeutic strategy for patients with Fabry disease, a lysosomal storage disease.

Immunogen information

Immunogen: galactosidase, alpha

Synonyms: Alpha galactosidase A, galactosidase, alpha

Calculated MW: 49kDa Uniprot ID: P06280

Application

Reactivity: Human

Tested Application: ELISA, IHC, WB, IF

Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>


1



Image:



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab00330(GLA Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00330(GLA antibody) at dilution of 1:1000