

anti- ALDOA antibody

Product Information

Catalog No.: FNab00301

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Plays a key role in glycolysis and gluconeogenesis. In addition, may also function as scaffolding protein(By similarity).

Immunogen information

Immunogen: aldolase A, fructose-bisphosphate

Synonyms: ALDOA; Aldolase A; Lung cancer antigen NY LU 1; Muscle type

aldolase

Observed MW: 39 kDa UniprotID: P04075

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB, IHC, IF, IP, FC

Recommended dilution: WB: 1:500-1:5000; IHC: 1:20-1:200; IP: 1:500-1:5000; IF: 1:10-1:100

Image:


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)


Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com

1





Immunohistochemistry of paraffin-embedded human skeletal muscle using FNab00301(ALDOA antibody) at dilution of 1:50




IP Result of anti-ALDOA (IP:FNab00301, 3ug; Detection:FNab00301 1:1000) with HepG2 cells lysate 6000ug.

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





HeLa cells were subjected to SDS PAGE followed by western blot with FNab00301(ALDOA antibody) at dilution of 1:500

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel: (0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com