

anti- ALG5 antibody

Product Information

| Catalog No.: | FNab00307 |
|---------------|---|
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | \geq 95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |
| Background | |

None

Immunogen information

| Immunogen: | asparagine-linked glycosylation 5, dolichyl-phosphate beta- glucosyltransferase homolog(S. cerevisiae) |
|--------------|---|
| Synonyms: | ALG5, DolP glucosyltransferase, bA421P11.2, RP11 421P11.2 |
| Observed MW: | 37 kDa |
| UniprotID : | Q9Y673 |

Application


| Reactivity: | Human, Mouse, Rat | | |
|---|-------------------|--|--|
| Tested Application: | ELISA, WB, IHC | | |
| Recommended dilution: WB: 1:200-1:2000; IHC: 1:20-1:200 | | | |
| Image: | | | |

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Immunohistochemistry of paraffin-embedded human placenta using FNab00307(ALG5 antibody) at dilution of 1:50



human kidney tissue were subjected to SDS PAGE followed by western blot with FNab00307(ALG5 antibody) at dilution of 1:500

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>