

# anti- Alpha-2-antiplasmin antibody

#### **Product Information**

Catalog No.: FNab00339

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

None

### **Immunogen information**

Immunogen: serpin peptidase inhibitor, clade F(alpha-2 antiplasmin, pigment

epithelium derived factor), member 2

Synonyms: A2AP, AAP, Alpha 2 antiplasmin, Alpha 2 AP, ALPHA 2 PI, Alpha 2

plasmin inhibitor, Alpha-2-antiplasmin, API, PLI, Serpin F2, SERPINF2

Observed MW: 54 kDa Uniprot ID: P08697

# **Application**

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:2000

Image:


#### 1

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





human liver tissue were subjected to SDS PAGE followed by western blot with FNab00339(SERPINF2 antibody) at dilution of 1:500

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>