

anti- alpha-SNAP antibody

Product Information

Catalog No.: FNab00343

Size: 100µg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

SNAPs(soluble NSF-attachment proteins) may play an important role in the process of docking and fusion of vesicles to their target membranes, according to the 'SNARE hypothesis'. NAPA(Alpha-SNAP) was found in a wide range of tissues and act synergistically in intra-Golgi transport. This antibody recognizes a 33kd protein in SDS page.

Immunogen information

Immunogen: N-ethylmaleimide-sensitive factor attachment protein, alpha

Synonyms: Alpha-SNAP, NAPA, SNAP alpha, SNAPA

Observed MW: 34 kDa Uniprot ID: P54920

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA
Recommended dilution: None

1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>