

anti- alpha-SNAP antibody

Product Information

Catalog No.:	FNab00344
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

SNAPs(soluble NSF-attachment proteins) may play an important role in the process of docking and fusion of vesicles to their target membranes, according to the 'SNARE hypothesis'. NAPA(Alpha-SNAP) was found in a wide range of tissues and act synergistically in intra-Golgi transport. This antibody recognizes a 33kd protein in SDS page.

Immunogen information

Immunogen:	N-ethylmaleimide-sensitive factor attachment protein, alpha
Synonyms:	Alpha-SNAP, NAPA, SNAP alpha, SNAPA
Observed MW:	34 kDa
UniprotID :	P54920

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:2000


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com

Image:



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab00344(alpha-SNAP antibody) at dilution of 1:500

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com