

anti- ALPK1 antibody

Product Information

Catalog No.: FNab00346

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Kinase that recognizes phosphorylation sites in which the surrounding peptides have an alphahelical conformation.

Immunogen information

Immunogen: alpha-kinase 1

Synonyms: 8430410J10Rik, alpha kinase 1, Alpha protein kinase 1, ALPK1,

Chromosome 4 kinase, FLJ22670, KIAA1527, LAK

Observed MW: 139 kDa UniprotID: Q96QP1

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:200-1:2000; IHC: 1:20-1:200

Image:


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)


Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com

1





Immunohistochemistry of paraffin-embedded human hepatocirrhosis using FNab00346(ALPK1 antibody) at dilution of 1:50



COLO 320 cells were subjected to SDS PAGE followed by western blot with FNab00346(ALPK1 antibody) at dilution of 1:500