

anti- ALS2 antibody

Product Information

Catalog No.:	FNab00350
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May act as a GTPase regulator. Controls survival and growth of spinal motoneurons(By similarity).

Immunogen information

Immunogen:	amyotrophic lateral sclerosis 2(juvenile)
Synonyms:	ALS2, ALS2CR6, Alsin, ALSJ, FLJ31851, IAHSP, KIAA1563, PLSJ
Observed MW:	185 kDa
UniprotID :	Q96Q42

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using FNab00350(ALS2 Antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab00350(ALS2 antibody) at dilution of 1:500