

anti- ALS2CR8 antibody

Product Information

Catalog No.:	FNab00351
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	amyotrophic lateral sclerosis 2(juvenile) chromosome region, candidate 8
Synonyms:	ALS2CR8, Calcium response factor, CARF, DKFZp667N246, FLJ12675, FLJ21579, NYD SP24
Observed MW:	81 kDa
Uniprot ID :	Q8N187

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:IHC:	1:20-1:200
Image:	

1


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab00351(ALS2CR8 Antibody) at dilution of 1:50