

anti- ALX1 antibody

Product Information

Catalog No.: FNab00352

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Sequence-specific DNA-binding transcription factor that binds palindromic sequences within promoters and may activate or repress the transcription of a subset of genes(PubMed:9753625, PubMed:8756334). Most probably regulates the expression of genes involved in the development of mesenchyme-derived craniofacial structures. Early on in development, it plays a role in forebrain mesenchyme survival(PubMed:20451171). May also induce epithelial to mesenchymal transition(EMT) through the expression of SNAI1(PubMed:23288509).

Immunogen information

Immunogen: ALX homeobox 1

Synonyms: ALX homeobox 1, ALX homeobox protein 1, ALX1, CART 1, CART 1,

Cartilage homeoprotein 1

Observed MW: 37 kDa
UniprotID: Q15699

Application

Reactivity: Human, Mouse

1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com



Tested Application: ELISA, WB

Recommended dilution: WB: 1:200-1:2000

Image:



mouse embryo tissue were subjected to SDS PAGE followed by western blot with FNab00352(ALX1 Antibody) at dilution of 1:600

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel: (0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com