

# anti- ALX3 antibody

## **Product Information**

Catalog No.: FNab00353

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

Transcriptional regulator with a possible role in patterning of mesoderm during development.

## **Immunogen information**

Immunogen: ALX homeobox 3

Synonyms: Homeobox protein aristaless-like 3,Proline-rich transcription factor

ALX3

Observed MW: 37 kDa, observed 70 kDa

UniprotID: O95076

# **Application**

Reactivity: Human, Mouse

Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:5000

Image:


#### 1

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel: (0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab00353(ALX3 Antibody) at dilution of 1:1000