

anti- AMIGO3 antibody

Product Information

Catalog No.:	FNab00361
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May mediate heterophilic cell-cell interaction. May contribute to signal transduction through its intracellular domain(By similarity).

Immunogen information

Immunogen:	adhesion molecule with Ig-like domain 3
Synonyms:	ALI3, Alivin 3, AMIGO 3, AMIGO3, Amphoterin induced protein 3, KIAA1851
Observed MW:	55 kDa
Uniprot ID :	Q86WK7

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab00361(AMIGO3 antibody) at dilution of 1:300