

anti- Aminoacylase 1 antibody

Product Information

Catalog No.:	FNab00362
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

This gene encodes a cytosolic, homodimeric, zinc-binding enzyme that catalyzes the hydrolysis of acylated L-amino acids to L-amino acids and an acyl group, and has been postulated to function in the catabolism and salvage of acylated amino acids. This gene is located on chromosome 3p21.1, a region reduced to homozygosity in small-cell lung cancer (SCLC), and its expression has been reported to be reduced or undetectable in SCLC cell lines and tumors. The amino acid sequence of human aminoacylase-1 is highly homologous to the porcine counterpart, and this enzyme is the first member of a new family of zinc-binding enzymes. Mutations in this gene cause aminoacylase-1 deficiency, a metabolic disorder characterized by central nervous system defects and increased urinary excretion of N-acetylated amino acids. Alternative splicing of this gene results in multiple transcript variants. Read-through transcription also exists between this gene and the upstream ABHD14A (abhydrolase domain containing 14A) gene, as represented in GeneID:100526760. A related pseudogene has been identified on chromosome 18.

Immunogen information

Immunogen:	aminoacylase 1
Synonyms:	ACY 1, ACY1, ACY1D, ACYLASE, Aminoacylase 1
Observed MW:	46 kDa
Uniprot ID :	Q03154

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel :(0086)027-87384275 F	ax: (0086)027-87800889
----------------------------	------------------------

www.fn-test.com



Application

Reactivity:Human, Mouse, RatTested Application:ELISA, WB, IHC, IFRecommended dilution:WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:50 - 1:200Image:



----70kd ----55kd ----40kd ----35kd ----25kd Immunohistochemistry of paraffin-embedded human esophageal cancer using FNab00362(ACY1 Antibody) at dilution of 1:50

K-562 cells were subjected to SDS PAGE followed by western blot with FNab00362(ACY1 Antibody) at dilution of 1:600

2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com