

anti- AMPD1-Specific antibody

Product Information

Catalog No.:	FNab00370	
Size:	100µg	
Form:	liquid	
Purification:	Immunogen affinity purified	
Purity:	\geq 95% as determined by SDS-PAGE	
Host:	Rabbit	
Clonality:	polyclonal	
Clone ID:	None	
IsoType:	IgG	
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)	

Background

AMPD1, also named as MAD and MADA, belongs to the adenosine and AMP deaminases family. AMP deaminase plays a critical role in energy metabolism. It catalyzes the reaction: AMP + H2O = IMP + NH3. Defects in AMPD1 are the cause of adenosine monophosphate deaminase deficiency muscle type(AMPDDM) which is a metabolic disorder resulting in exercise-related myopathy. The antibody is specific to AMPD1.

Immunogen information

Immunogen:	adenosine monophosphate deaminase 1(isoform M)
Synonyms:	AMPD1, MAD, Myoadenylate deaminase, AMP deaminase isoform M, AMP deaminase 1
Observed MW:	80 kDa
UniprotID :	P23109

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com

1



Recommended dilution:WB: 1:500-1:5000 Image:



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab00370(AMPD1-Specific antibody) at dilution of 1:300

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>