

anti- ANGPTL2 antibody

Product Information

Catalog No.:	FNab00397
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Induces sprouting in endothelial cells through an autocrine and paracrine action.

Immunogen information

Immunogen:	angiopoietin-like 2
Synonyms:	Angiopoietin like 2, Angiopoietin like protein 2, Angiopoietin related protein 2, ANGPTL2, ARP2, HARP
Observed MW:	57 kDa
Uniprot ID :	Q9UKU9

Application

Reactivity:	Human, Mouse, Rat	
Tested Application:	ELISA, WB, IHC, IP	
Recommended dilution: WB: 1:500-1:5000; IHC: 1:200-1:2000; IP: 1:20-1:200		
Image:		


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889


www.fn-test.com





Immunohistochemistry of paraffin-embedded human heart using FNab00397(ANGPTL2 antibody) at dilution of 1:50

IP Result of anti-ANGPTL2 (IP:FNab00397, 4ug; Detection:FNab00397 1:500) with mouse testis tissue lysate 4800ug.




2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com





3

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)