

# anti- ANGPTL7 antibody

#### **Product Information**

Catalog No.: FNab00401

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

None

# Immunogen information

Immunogen: angiopoietin-like 7

Synonyms: Angiopoietin like 7, Angiopoietin like factor, Angiopoietin like protein 7,

Angiopoietin related protein 7, ANGPTL7, AngX, CDT6, dJ647M16.1,

RP4 647M16.2

Observed MW: 45 kDa Uniprot ID: 043827

## **Application**

Reactivity: Human, Mouse

Tested Application: ELISA, WB

Recommended dilution: WB: 1:200-1:2000

Image:


### 1

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





MCF7 cells were subjected to SDS PAGE followed by western blot with FNab00401(ANGPTL7 antibody) at dilution of 1:800