

anti- ANKRD11 antibody

Product Information

Catalog No.: FNab00405

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

ANKRD11, also named as ANCO1, may recruit HDACs to the p160 coactivators/nuclear receptor complex to inhibit ligand-dependent transactivation. This antibody is specific to ANKRD11, ANKRD11 exists 150-170 kDa isoform.

Immunogen information

Immunogen: ankyrin repeat domain 11

Synonyms: ANCO 1, ANCO1, ANKRD11, ankyrin repeat domain 11, LZ16, T13

Observed MW: 298 kDa, Observed 160 kDa

UniprotID: Q6UB99

Application

Reactivity: Human

Tested Application: ELISA, WB

Recommended dilution: WB: 1:200-1:2000

Image:


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com

1





BxPC-3 cells were subjected to SDS PAGE followed by western blot with FNab00405(ANKRD11 antibody) at dilution of 1:500

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel: (0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com