

anti- ANKRD53 antibody

Product Information

Catalog No.:	FNab00416
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Required for normal progression through mitosis. Involved in chromosome alignment and cytokinesis via regulation of microtubules polymerization.

Immunogen information

Immunogen:	ankyrin repeat domain 53
Synonyms:	ANKRD53, ankyrin repeat domain 53, FLJ12056, FLJ36160
Observed MW:	60 kDa
UniprotID :	Q8N9V6

Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000; IF: 1:50-1:500; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab00416(ANKRD53 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00416(ANKRD53 Antibody) at dilution of 1:1000