

anti- ANKRD55 antibody

Product Information

Catalog No.:	FNab00418
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	ankyrin repeat domain 55
Synonyms:	ANKRD55, ankyrin repeat domain 55, FLJ11795
Observed MW:	68 kDa
UniprotID :	Q3KP44

Application


Reactivity:	Human
Tested Application:	ELISA, IHC, IF, WB
Recommended dilution:	WB: 1:200-1:2000; IF: 1:20-1:200; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human cervical cancer tissue slide using FNab00418(ANKRD55 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00418(ANKRD55 Antibody) at dilution of 1:600