

anti- ANKRD57 antibody

Product Information

Catalog No.:	FNab00419
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	ankyrin repeat domain 57
Synonyms:	ANKRD57, ankyrin repeat domain 57, C2orf26
Observed MW:	56 kDa, Observed 70 kDa
UniprotID :	Q53LP3

Application


Reactivity:	Human
Tested Application:	ELISA, IP, IHC
Recommended dilution:	IP:1:500-1:5000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human cerebellum tissue slide using FNab00419(ANKRD57 Antibody) at dilution of 1:50



IP Result of anti-ANKRD57 (IP:FNab00419, 3ug; Detection:FNab00419 1:1000) with HeLa cells lysate 3000ug.