

anti- ANKRD7 antibody

Product Information

Catalog No.: FNab00421

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen: ankyrin repeat domain 7

Synonyms: ANKRD7, ankyrin repeat domain 7, Testis specific protein TSA806,

TSA806

Observed MW: 29 kDa UniprotID: Q92527

Application

Reactivity: Human

Tested Application: ELISA, WB

Recommended dilution: WB: 1:200-1:1000

Image:


1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





human testis tissue were subjected to SDS PAGE followed by western blot with FNab00421(ANKRD7 Antibody) at dilution of 1:300

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com