

anti- ANKS1B antibody

Product Information

Catalog No.: FNab00422

Size: 100µg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Isoform 2 may participate in the regulation of nucleoplasmiccoilin protein interactions in neuronal and transformed cells. Isoform 3 can regulate global protein synthesis by altering nucleolar numbers. Isoform 4 may play a role as a modulator of APP processing. Overexpression can down-regulate APP processing.

Immunogen information

Immunogen: ankyrin repeat and sterile alpha motif domain containing 1B

Synonyms: AIDA, AIDA 1, ANKS1B, ANKS2, cajalin 2, E2A PBX1 associated

protein, EB 1, EB1

Observed MW: 138 kDa, Observed 70 kDa

UniprotID: Q7Z6G8

Application

Reactivity: Human, Mouse

Tested Application: ELISA, WB, IF, IHC

Recommended dilution: WB: 1:500-1:5000; IP: 1:500-1:5000; IHC: 1:20-1:200; IF: 1:20-1:200

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com


1



Image:



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab00422(ANKS1B Antibody) at dilution of 1:50



A549 cells were subjected to SDS PAGE followed by western blot with FNab00422(ANKS1B Antibody) at dilution of 1:1000