

anti- ANKS3 antibody

Product Information

Catalog No.: FNab00423

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen: ankyrin repeat and sterile alpha motif domain containing 3

Synonyms: ANKS3, FLJ32345, FLJ32767, KIAA1977

Observed MW: 72 kDa, 64 kDa

UniprotID: Q6ZW76

Application

Reactivity: Human

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500-1:5000; IHC: 1:20-1:200

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab00423(ANKS3 Antibody) at dilution of 1:50



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab00423(ANKS3 Antibody) at dilution of 1:1000