

anti- ANKZF1 antibody

Product Information

Catalog No.: FNab00424

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Involved in the endoplasmic reticulum(ER)-associated degradation(ERAD) pathway.

Immunogen information

Immunogen: ankyrin repeat and zinc finger domain containing 1

Synonyms: ANKZF1, FLJ10415, FLJ13144, Zinc finger protein 744, ZNF744

Observed MW: 81 kDa
UniprotID: Q9H8Y5

Application

Reactivity: Human

Tested Application: ELISA, WB

Recommended dilution: WB: 1:1000-1:10000

Image:


1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





BxPC-3 cells were subjected to SDS PAGE followed by western blot with FNab00424(ANKZF1 Antibody) at dilution of 1:2000