

# anti-Annexin IV antibody

#### **Product Information**

Catalog No.: FNab00434

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

### **Background**

Annexin IV /placental anticoagulant protein II(ANX4) is a member of the annexin or lipocortin family of calcium-dependent phospholipid binding proteins, and is located at 2p13. Although their functions are still not clearly defined, several members of the annexin family have been implicated in membrane-related events along exocytotic and endocytotic pathways. ANX4 shares 45 to 59% identity with other members and possibly interacts with ATP. ANX4 has in vitro anticoagulant activity and also inhibits phospholipase A2 activity. ANX4 is almost exclusively expressed in epithelial cells. In western blotting, we get 36 kDa, 39 kDa(with modification), and 60 kDa(an unspecifc band).

#### **Immunogen information**

Immunogen: annexin A4

Synonyms: 35 beta calcimedin, Annexin 4, annexin A4, Annexin IV, ANX4,

ANXA4, Chromobindin 4, DKFZp686H02120, Endonexin I, Lipocortin

IV, P32.5, PAP II, PIG28, PP4 X, Protein II, ZAP36

Observed MW: 36 kDa UniprotID: P09525

#### 1

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>




## **Application**


Reactivity: Human, Mouse, Rat Tested Application: ELISA, IHC, WB

Recommended dilution: WB: 1:500-1:5000; IHC: 1:20-1:200

Image:



Immunohistochemistry of paraffin-embedded human stomach tissue slide using FNab00434(Annexin IV Antibody) at dilution of 1:50



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with FNab00434(Annexin IV Antibody) at dilution of 1:1000