

anti- APH1A antibody

Product Information

Catalog No.:	FNab00481	
Size:	100µg	
Form:	liquid	
Purification:	Immunogen affinity purified	
Purity:	\geq 95% as determined by SDS-PAGE	
Host:	Rabbit	
Clonality:	polyclonal	
Clone ID:	None	
IsoType:	IgG	
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)	

Background

Essential subunit of the gamma-secretase complex, an endoprotease complex that catalyzes the intramembrane cleavage of integral proteins such as Notch receptors and APP(beta-amyloid precursor protein). It probably represents a stabilizing cofactor for the presenilin homodimer that promotes the formation of a stable complex.

Immunogen information

Immunogen:	anterior pharynx defective 1 homolog A
Synonyms:	APH 1A, Aph 1alpha, APH1A, CGI 78, Gamma secretase subunit APH 1A, PSF, UNQ579/PRO1141
Observed MW:	30 kDa
Uniprot ID :	Q96BI3

Application

Reactivity:	Human, Mouse, Rat		
Tested Application:	ELISA, WB		
Recommended dilution:WB: 1:500 - 1:1000			


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206) Tel :(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com

1



Image:



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00481(APH1A antibody) at dilution of 1:500

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com