

## anti- APOA4 antibody

### Product Information

Catalog No.:	FNab00490
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	3E8
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

May have a role in chylomicrons and VLDL secretion and catabolism. Required for efficient activation of lipoprotein lipase by ApoC-II; potent activator of LCAT. ApoA-IV is a major component of HDL and chylomicrons.

### Immunogen information

Immunogen:	apolipoprotein A-IV
Synonyms:	None
Observed MW:	45-46 kDa
UniprotID :	P06727

### Application


Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:1000-1:5000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



L02 cells were subjected to SDS PAGE followed by western blot with FNab00490( APOA4 Antibody) at dilution of 1:2000