

# anti- APOBEC3B antibody

## **Product Information**

Catalog No.:	FNab00493	
Size:	100µg	
Form:	liquid	
Purification:	Immunogen affinity purified	
Purity:	$\geq$ 95% as determined by SDS-PAGE	
Host:	Rabbit	
Clonality:	polyclonal	
Clone ID:	None	
IsoType:	IgG	
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)	

### Background

APOBEC3B has 3 isoforms produced by alternative splicing with the MW of 46, 29, 43 kDa. It can also exist as a dimer and can interact with APOBEC3G.

#### **Immunogen information**

Immunogen:	apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3B
Synonyms:	None
Observed MW:	46 kDa
UniprotID :	Q9UH17

# Application


Reactivity:	Human		
Tested Application:	ELISA, WB, IHC		
Recommended dilution: WB: 1:200-1:1000; IHC: 1:20-1:200			
Image:			

#### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab00493(APOBEC3B Antibody) at dilution of 1:200 heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



Hela cells were subjected to SDS PAGE followed by western blot with FNab00493(APOBEC3B Antibody) at dilution of 1:300

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>