

anti- APRT antibody

Product Information

Catalog No.:	FNab00518
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Adenine phosphoribosyltransferase belongs to the purine/pyrimidine phosphoribosyltransferase family. A conserved feature of this gene is the distribution of CpG dinucleotides. This enzyme catalyzes the formation of AMP and inorganic pyrophosphate from adenine and 5-phosphoribosyl-1-pyrophosphate (PRPP). It also produces adenine as a by-product of the polyamine biosynthesis pathway. A homozygous deficiency in this enzyme causes 2, 8-dihydroxyadenine urolithiasis. Two transcript variants encoding different isoforms have been found for this gene.

Immunogen information

Immunogen:	adenine phosphoribosyltransferase
Synonyms:	None
Observed MW:	21 kDa
Uniprot ID :	P07741

Application

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889


www.fn-test.com



Tested Application: ELISA, WB, IHC Recommended dilution: WB: 1:500 - 1:2000; IHC: 1:50 - 1:200 Image:



Immunohistochemistry of paraffin-embedded mouse lung using FNab00518(APRT antibody) at dilution of 1:50



human blood tissue were subjected to SDS PAGE followed by western blot with FNab00518(APRT antibody) at dilution of 1:1000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com