

anti- ARHGAP21 antibody

Product Information

| Catalog No.: | FNab00550 |
|---------------|---|
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | \geq 95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Background

Functions as a GTPase-activating protein(GAP) for RHOA and CDC42. Downstream partner of ARF1 which may control Golgi apparatus structure and function. Also required for CTNNA1 recruitment to adherens junctions.

Immunogen information

| Immunogen: | Rho GTPase activating protein 21 |
|--------------|----------------------------------|
| Synonyms: | ARHGAP10, KIAA1424 |
| Observed MW: | 250 kDa |
| UniprotID : | Q5T5U3 |


Application

| Reactivity: | Human, Rat | | |
|---|--------------------|--|--|
| Tested Application: | ELISA, WB, IF, IHC | | |
| Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200 | | | |
| Image: | | | |


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan,Hubei, China(430206)Tel :(0086)027-87384275Fax: (0086)027-87800889 www.fn-test.com





Immunohistochemistry of paraffin-embedded human ovary tissue slide using FNab00550(ARHGAP21 Antibody) at dilution of 1:50



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab00550(ARHGAP21 Antibody) at dilution of 1:600

Wuhan Fine Biotech Co., Ltd. B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>