

# anti- ARHGAP25 antibody

## **Product Information**

Catalog No.:	FNab00551
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

# Background

ARHGAPs, such as ARHGAP25, encode negative regulators of Rho GTPases (see ARHA; MIM 165390), which are implicated in actin remodeling, cell polarity, and cell migration (Katoh and Katoh, 2004 [PubMed 15254788]).

### Immunogen information

Immunogen:	Rho GTPase activating protein 25
Synonyms:	KIAA0053
Observed MW:	70 kDa
Uniprot ID :	P42331

# Application


Reactivity:	Human		
Tested Application:	ELISA, WB, IHC		
Recommended dilution: WB: 1:500 - 1:2000; IHC: 1:50 - 1:200			
Image:			

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206) Tel :( 0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com

1





Immunohistochemistry of paraffin-embedded human gastric using FNab00551( ARHGAP25 antibody) at dilution of 1:200



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab00551(ARHGAP25 antibody) at dilution of 1:1000

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com