

anti- ARHGEF18 antibody

Product Information

Catalog No.:	FNab00561
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Acts as guanine nucleotide exchange factor(GEF) for RhoAGTPases. Its activation induces formation of actin stress fibers. Also act as a GEF for RAC1, inducing production of reactive oxygen species(ROS). Does not act as a GEF for CDC42. The G protein beta-gamma(Gbetagamma) subunits of heterotrimeric G proteins act as activators, explaining the integrated effects of LPA and other G-protein coupled receptor agonists on actin stress fiber formation, cell shape change and ROS production. Required for EPB41L4B-mediated regulation of the circumferential actomyosin belt in epithelial cells(PubMed:22006950).

Immunogen information

Immunogen:	rho/rac guanine nucleotide exchange factor(GEF) 18
Synonyms:	KIAA0521
Observed MW:	120-130 kDa
UniprotID :	Q6ZSZ5

Application

Reactivity: Human, Mouse


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>



Tested Application: ELISA, WB, IHC Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200 Image:



-55kd

Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab00561(ARHGEF18 antibody) at dilution of 1:50

mouse testis tissue were subjected to SDS PAGE -170kd 1:300 -130kd -100kd -70kd

followed by western blot with FNab00561(ARHGEF18 antibody) at dilution of

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com