

anti- ARL8A antibody

Product Information

Catalog No.: FNab00581

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24

months (Avoid repeated freeze / thaw cycles.)

Background

May play a role in lysosomes motility. Alternatively, may play a role in chromosome segregation(By similarity).

Immunogen information

Immunogen: ADP-ribosylation factor-like 8A

Synonyms: ARL10B, ARL8A, GIE2

Observed MW: 22 kDa Uniprot ID: Q96BM9

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB, IP

Recommended dilution: WB: 1:500-1:1000; IP: 1:500-1:1000

Image:


1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab00581(ARL8A Antibody) at dilution of 1:600