

anti- ARL8B antibody

Product Information

Catalog No.: FNab00582

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24

months (Avoid repeated freeze / thaw cycles.)

Background

May play a role in lysosome motility(PubMed:16537643, PubMed:25898167). May play a role in chromosome segregation(PubMed:15331635).

Immunogen information

Immunogen: ADP-ribosylation factor-like 8B

Synonyms: ARL10C, ARL8B, Gie1

Observed MW: 22 kDa Uniprot ID: Q9NVJ2

Application

Reactivity: Human, Mouse

Tested Application: ELISA, IHC, WB

Recommended dilution: WB: 1:500-1:1000; IHC: 1:20-1:200

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human stomach tissue slide using FNab00582(ARL8B Antibody) at dilution of 1:50



fetal human brain tissue were subjected to SDS PAGE followed by western blot with FNab00582(ARL8B Antibody) at dilution of 1:600